



ALTA750G[®]

X-RAY TUBE HOUSING

PRODUCT DESCRIPTION

The ALTA750G[®] is an x-ray tube housing assembly specifically designed for use with CT Scanners. This tube housing assembly consists of the ALTA750G[®] Tube loaded into the Richardson ALTA750G[®] Housing or reloaded into the Varex* B-805H Housing. The Richardson G Heat Exchanger or the Varex* HE-978 Exchanger can be used with this tube housing assembly.

INTENDED USE

The ALTA750G[®] X-Ray Tube Housing Assemblies are designed to emit ionizing radiation and are intended to be used as a component of a CT system which is used for diagnostic and interventional X-Ray applications.

INCLUDED INFORMATION AND SPECIFICATIONS



- Tube Specifications
- Housing Assembly Specifications
- Volumetric / Helical Scan Ratings
- Cathode Emission Characteristics
- Housing Diagram
- Housing Wiring
- Disposal Information

Originally written in English.


Additional copies and alternate language versions available upon request from techdata@rell.com

*This product is not affiliated with, endorsed by, or sponsored by Varex Imaging.

TUBE SPECIFICATIONS

| | | |
|--|-----------|-----------|
| Nominal X-ray Tube Voltage | kV | 140 |
| Anode Diameter | mm | 200 |
| Anode Material | | ReW-TZM-C |
| Anode Angle | Degrees | 7 |
| Nominal Focal Spot – Small  IEC 60336 Loading Factor 120kV x 200mA | IEC 60336 | 0.9 x 0.8 |
| Nominal Focal Spot – Large  IEC 60336 Loading Factor 120kV x 200mA | IEC 60336 | 1.6 x 1.4 |
| Anode Heat Content (Maximum) | MJ | 5.4 |
| Nominal Anode Input Power – Large IEC 60613 | kW | 72 |
| Nominal Anode Input Power – Small IEC 06013 | kW | 42 |
| Nominal CT Scan Power Index – Large IEC 60613 | kW | 69.5 |
| Nominal CT Scan Power Index – Small IEC 60613 | kW | 42 |
| Anode Heat Dissipation (Maximum) | W | 12,000 |
| Maximum Filament Current - Large | A | 5 |
| Maximum Filament Voltage - Large | V | 14.4 |
| Maximum Filament Current - Small | A | 4.8 |
| Maximum Filament Voltage - Small | V | 12.3 |

HOUSING ASSEMBLY SPECIFICATIONS

| | | |
|--|-----------|-----------|
| Maximum Heat Content | MJ | 3.6 |
| Maximum Continuous Heat Dissipation | kW | 4.0 |
| Maximum Housing Temperature | Degrees C | 78 |
| Permanent Filtration  IEC 60522 | mm AL | 1.0 |
| Temperature Limits for Transport and Storage | Degrees C | -20 to 75 |
| Temperature Limits for Operation | Degrees C | 5 to 40 |
| Weight of Assembly | kg | 60.5 |
| Leakage Radiation @140 kV, 29 mA, 1 m | mGy/hr | 0.57 |

ADDITIONAL HOUSING ASSEMBLY SPECIFICATIONS

| |
|--|
| Humidity Limits for Transport and Storage: 10% to 90% RH |
| Pressure Limits for Transport and Storage: 70 to 106 kPa |
| Humidity Limits for Normal Operation: 40 to 80% RH |
| Pressure Limits for Normal Operation: 70 to 106 kPa |
| Classification per IEC 60601-1: Class 1 Type B |
| Degree of Protection Against Ingress of Water is IPX0 |
| Mode of Operation: Intermittent (non-continuous) |
| Device Classification: U.S FDA = Class 1, EU = Class IIb |

VOLUMETRIC / HELICAL SCAN RATINGS IEC 60613

\varnothing 50 Hz
 0.9 x 0.8 Focal Spot
 7 Degree Target

| Volume Scan Time (Seconds) | Maximum Allowed Tube Current (mA) as a Function of the Following Starting Heat Storage and Tube Voltages | | | | | | | | |
|----------------------------|---|--------|--------|---------------------------|--------|--------|---------------------------|--------|--------|
| | Starting Heat Storage 40% | | | Starting Heat Storage 55% | | | Starting Heat Storage 70% | | |
| | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV |
| 4 | 300 | 250 | 225 | 300 | 250 | 225 | 300 | 250 | 225 |
| 10 | 300 | 250 | 225 | 300 | 250 | 225 | 300 | 250 | 225 |
| 15 | 300 | 250 | 225 | 300 | 250 | 225 | 300 | 250 | 225 |
| 20 | 300 | 250 | 225 | 300 | 250 | 225 | 300 | 250 | 225 |
| 30 | 300 | 250 | 225 | 300 | 250 | 225 | 300 | 250 | 225 |
| 45 | 300 | 250 | 225 | 300 | 250 | 225 | 300 | 250 | 225 |
| 60 | 300 | 250 | 225 | 300 | 250 | 225 | 250 | 200 | 175 |
| 75 | 300 | 250 | 225 | 300 | 250 | 225 | 225 | 175 | 150 |
| 80 | 300 | 250 | 225 | 300 | 250 | 225 | 200 | 175 | 150 |
| 90 | 300 | 250 | 225 | 275 | 225 | 200 | 200 | 150 | 150 |

\varnothing 50 Hz
 1.6 x 1.4 Focal Spot
 7 Degree Target

| Volume Scan Time (Seconds) | Maximum Allowed Tube Current (mA) as a Function of the Following Starting Heat Storage and Tube Voltages | | | | | | | | |
|----------------------------|---|--------|--------|---------------------------|--------|--------|---------------------------|--------|--------|
| | Starting Heat Storage 40% | | | Starting Heat Storage 55% | | | Starting Heat Storage 70% | | |
| | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV |
| 4 | 670 | 560 | 500 | 670 | 560 | 500 | 670 | 560 | 500 |
| 10 | 670 | 560 | 500 | 670 | 560 | 500 | 670 | 560 | 490 |
| 15 | 670 | 560 | 500 | 670 | 560 | 500 | 640 | 530 | 470 |
| 20 | 670 | 560 | 500 | 670 | 560 | 500 | 610 | 510 | 450 |
| 30 | 600 | 500 | 440 | 600 | 500 | 440 | 440 | 360 | 320 |
| 45 | 540 | 450 | 400 | 480 | 400 | 350 | 320 | 270 | 240 |
| 60 | 450 | 370 | 330 | 380 | 310 | 280 | 260 | 220 | 190 |
| 75 | 410 | 340 | 300 | 310 | 260 | 230 | 230 | 190 | 170 |
| 80 | 380 | 320 | 280 | 300 | 250 | 220 | 220 | 180 | 160 |
| 90 | 350 | 290 | 260 | 270 | 230 | 200 | 200 | 170 | 150 |



\varnothing 100 Hz
 0.9 x 0.8 Focal Spot
 7 Degree Target

| Volume Scan Time (Seconds) | Maximum Allowed Tube Current (mA) as a Function of the Following Starting Heat Storage and Tube Voltages | | | | | | | | |
|----------------------------|---|--------|--------|---------------------------|--------|--------|---------------------------|--------|--------|
| | Starting Heat Storage 40% | | | Starting Heat Storage 55% | | | Starting Heat Storage 70% | | |
| | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV |
| 4 | 425 | 350 | 300 | 425 | 350 | 300 | 425 | 350 | 300 |
| 10 | 425 | 350 | 300 | 425 | 350 | 300 | 425 | 350 | 300 |
| 15 | 425 | 350 | 300 | 425 | 350 | 300 | 425 | 350 | 300 |
| 20 | 425 | 350 | 300 | 425 | 350 | 300 | 425 | 350 | 300 |
| 30 | 425 | 350 | 300 | 425 | 350 | 300 | 400 | 325 | 300 |
| 45 | 425 | 350 | 300 | 425 | 350 | 300 | 300 | 250 | 225 |
| 60 | 425 | 350 | 300 | 375 | 300 | 275 | 250 | 200 | 175 |
| 75 | 400 | 325 | 300 | 300 | 250 | 225 | 225 | 175 | 150 |
| 80 | 375 | 300 | 275 | 300 | 250 | 225 | 200 | 175 | 150 |
| 90 | 350 | 275 | 250 | 275 | 225 | 200 | 200 | 150 | 150 |


\varnothing 100 Hz
 1.6 x 1.4 Focal Spot
 7 Degree Target

| Volume Scan Time (Seconds) | Maximum Allowed Tube Current (mA) as a Function of the Following Starting Heat Storage and Tube Voltages | | | | | | | | |
|----------------------------|---|--------|--------|---------------------------|--------|--------|---------------------------|--------|--------|
| | Starting Heat Storage 40% | | | Starting Heat Storage 55% | | | Starting Heat Storage 70% | | |
| | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV | 100 kV | 120 kV | 135 kV |
| 4 | 720 | 600 | 530 | 720 | 600 | 530 | 720 | 600 | 530 |
| 10 | 720 | 600 | 530 | 720 | 600 | 530 | 720 | 600 | 530 |
| 15 | 720 | 600 | 530 | 720 | 600 | 530 | 720 | 600 | 530 |
| 20 | 720 | 600 | 530 | 720 | 600 | 530 | 610 | 510 | 450 |
| 30 | 600 | 500 | 440 | 600 | 500 | 440 | 440 | 360 | 320 |
| 45 | 540 | 450 | 400 | 480 | 400 | 350 | 320 | 270 | 240 |
| 60 | 450 | 370 | 330 | 380 | 310 | 280 | 260 | 220 | 190 |
| 75 | 410 | 340 | 300 | 310 | 260 | 230 | 230 | 190 | 170 |
| 80 | 380 | 320 | 280 | 300 | 250 | 220 | 220 | 180 | 160 |
| 90 | 350 | 290 | 260 | 270 | 230 | 200 | 200 | 170 | 150 |


CATHODE EMISSION CHARACTERISTICS IEC 60613




ANODE HEATING AND COOLING CURVES IEC 60613




TUBE HOUSING ASSEMBLY HEATING AND COOLING CURVES



HOUSING DIAGRAM



HOUSING WIRING



DISPOSAL INFORMATION

Take back, proper disposal and recovery of medical devices takes place in accordance with European WEEE directive and the requirements of national legislation.

The x-ray tube contains beryllium. The x-ray tube housing assembly contains lead for radiation shielding and mineral oil. The x-ray tube and x-ray tube housing must not be disposed in domestic or industrial waste; they must be disposed in accordance with local regulation.

The tube and housing assembly may be returned to Richardson Electronics, Ltd. for proper disposal.

Richardson Electronics, Ltd strives to be environmentally conscious. Select materials and components are recycled. Controls are in place to assure that all product meet specifications and safety requirements.

Richardson Electronics, Ltd. | 40W267 Keslinger Road P.O. Box 393 | LaFox, IL 60147-0393 | (630) 208-2200